

Shoulder Suture Passer

- ✓ Designed for rotator cuff and labral repair
- ✓ Aggressive needle angle for improved suture passing
- ✓ Engineered to fire only when intended
- ✓ Ergonomic, lightweight with engageable and disengageable ratchet
- ✓ Options to load the suture thread through the middle or from the side
- ✓ Single use nitinol needle providing for super-elasticity and strength
- ✓ Sterile packed and non-sterile packed needle options
- ✓ US FDA approved



Insert needle through the luer lock



Lock needle into the suture passer





Put suture in place



Grasp tissue



Pass and retrieve suture



Suture Passers

Item Number	Opening	Handle
STS-ART-400-170SP	Middle/Straight	Without Ratchet
STS-ART-401-170SP	Middle/Straight	With Ratchet
STS-ART-402-170SP	Sideward	Without Ratchet
STS-ART-403-170SP	Sideward	With Ratchet



Middle/Straight



Sideward



Single Use Nitinol Needles

Item Number	Description
STS-ART-430-170NN	Non-sterile packed
STS-ART-431-170NN*	Sterile packed, 6 per Box

* Coming soon